

Disassembly



STEP 1

Disassemble

Before cleaning, disassemble all three parts of the Provale Cup for proper cleaning. Grasp the thumb tab on the lid to remove the lid and inner cup assembly from the outer cup. Hold the inner cup "cone" in one hand and grasp the thumb tab on the lid to separate the two parts.

STEP 2

Wash

Wash each cup part in hot soapy water and rinse thoroughly or place each part on the top rack of a commercial or residential dishwasher. Parts are dishwasher safe up to 180 degrees Fahrenheit.

STEP 3

Air Dry

Place the parts on a clean surface or drying rack and allow them to air dry completely before reassembly.

STEP 4

Reassemble

Once dry, carefully reassemble your Provale Cup. For help with this, refer to the assembly section below.

STEP 5

Check Functionality

Make sure the inner cup is seated properly within the outer cup and the lid fits snugly. All vent holes must be clear and unobstructed for the cup to function properly.

Assembly



STEP 1

Attaching Handles

Determine if you would like to use handles. The cup can be used with one handle, both handles, or no handles. Slide each handle into one of the grooves on the side of the outer cup. Firmly pull the handle upward until a light snap is heard and the "ledge" of the handle is over the lip of the cup. Test the handles. If they are easily removed, repeat this step.

Once properly attached, handles cannot be removed. Attempting to remove handles will damage the cup and cause cup leakage.

STEP 2

Filling the Outer Cup

Once all Provale Cup parts are clean and dry, determine the correct cup size: the 5cc Modern Blue Cup or the 10cc Light Grey Cup. Confirm that the conical inner cup and lid that fits on top are the same color before proceeding. Fill the transparent cup base with thin liquid up to or below the fill line. The fill line is located on the inside of the outer cup (approximately 8 oz). Do not overfill, as this may cause leakage.

STEP 3

Inserting the Lid

Assemble the inner cup and the lid by snapping them together — an audible click should be heard. These parts can only be assembled one way.

STEP 4

Final Check

Place the conical inner cup and lid assembly into the transparent outer cup and press down on the edges. A slight audible click should be heard. The cup is now ready to use.

STEP 5

Refilling

To refill, remove the lid and inner cup assembly, then pour the liquid into the transparent base up to or below the 8 ounce fill line. Snap the lid and inner cup assembly back onto the transparent base.